- F. Manual RBC counts Example only
 - 1. Manual RBC counts are not done on EDTA whole blood, **HOWEVER**, cell counts on bloody body fluids (spinal fluid) may require making a RBC dilution of 1:200 with isotonic saline (which <u>prevents</u> lysis of RBCs). Red cells in the RBC counting area of 0.4mm² are counted (5 small squares in the center square of the hemocytometer on each side).
 - 2. Procedure
 - a. MIX EDTA whole blood specimen. Make a <u>1:200</u> dilution with saline and mix.
 - b. Plate dilution on each side of chamber and allow cells to settle 3 mins.
 - c. Count **4 outer squares** of <u>center sqmm</u> and the <u>center square of <u>center sqmm</u></u> (shaded squares on chamber) on <u>each side</u> of the hemocytometer using 40x objective and low light.
 - d. Total area counted = $0.4 \text{ mm}^2 \rightarrow \text{RBC}$ counting area.

					10.000	
	1111		1111	нш.	++	
1			1110			
	1111				1.1	1
	1111			****	11	
-				- -		
_						- F
					-	1
					-	
_						
		ШЦ		1111	+	_
					1.	
		III H				
					T	
					+ +	++

Center sqmm is divided into 25 squares; each square in center square is 1/25mm²;

5/25mm² counted on <u>each</u> side = 10/25mm² OR 1/5mm² + 1/5mm² = 2/5mm² = 0.4mm².

e. Calculations:

 $\frac{\text{RBC/cmm} = \frac{\text{\#cells (both sides) x 200}}{0.1 \text{ mm x } 0.4 \text{ sqmm}} \text{ or } \frac{\text{total \#cells x 200 x 10 mm}}{0.4 \text{ sqmm}}$

Report to nearest hundreth.

f.	Example:		
	Side 1 305	$620 \ge 200 \ge 10 = 3.10$ m	illion/cmm or $3.10 \ge 10^6/\text{uL}$
	Side 2 315	0.4	or 3.10 M/uL
			or 3.10×10^{12} /L (SI units)

SEE Calculations sheet for formulas, calculation examples and reporting units.

Complete Blood Count with Differential			Equivalent, Conventional (preferred) & SI Units
CBC paramete WBC RBC HGB HCT MCV MCH MCHC RDW PLT	8.8 5.01 15.0 44.9 86.6 29.0 33.5 11.5 211	thousands/uL millions/cmm g/dl % fl pg % % thousands/cmm	Equivalent units •cmm (cubic millimeters) = mm ³ = uL (microliters) = mcL (conventional units) •thousands/uL = $x10^{3}$ /cmm = $x10E3$ /uL = K/cmm <i>Examples:</i> WBC 6,400/uL = 6.4×10^{3} /cmm = $6.4 \times 10E3$ /uL = 6.4 K/cmm PLT 250,000/uL = 250×10^{3} /cmm = $250 \times 10E3$ /uL = 250 K/cmm •millions/uL = $x10^{6}$ /cmm = $x10E6$ /uL = M/cmm <i>Example:</i> RBC 3.45 million/uL = 3.45×10^{6} /cmm = $3.45 \times 10E6$ /uL = 3.45 M/cmm
Differential (W Neutrophils Lymphocytes Monocytes Eosinophils Basophils	BC types 56 34 8 1 1) % % % %	To convert conventional units (cmm or uL) to SI units (L), use a factor of x10⁶ <i>Examples:</i> WBC 6.4×10^{3} /uL (x10⁶) = 6.4×10^{9} /L (SI units) RBC 3.45×10^{6} /uL (x10⁶) = RBC 3.45×10^{12} /L (SI) PLT 250×10^{3} /uL (x10⁶) = 250×10^{9} /L (SI = System of International)

CLS416 Clinical Hematology